

Sofclear

COMFORT

with *BioMoist*

Try *BioMoist* today and
feel the difference !

By combining the benefits of Radial Edge technology with the exclusive properties of BioMoist, Gelflex has designed a contact lens that provides superior all day comfort.

The BioMoist formula is a hydrating solution containing hydroxypropyl methylcellulose (HPMC) and polyethylene glycol (PEG), which together mimics the natural tear film of the eye. The combination of these superior wetting agents provides a soft, silky feel which helps to optimise comfort upon initial insertion of the lens.

The specially formulated BioMoist solution releases moisture when it is needed most, resulting in superior comfort throughout the entire day.

BioMoist also contains unique cleaning agents that help to resist the build-up of protein and lipids on the lens surface, further enhancing the overall wearing experience.

The addition of BioMoist to our comfort range of Sofclear lenses demonstrates our commitment to producing an exceptionally comfortable contact lens.



ONE DAY LENSES
30 Pack

Base curve	8.60
Diameter	14.30
Powers available	-0.50 to -6.00 in 0.25 steps -6.50 to -10.00 in 0.50 steps +0.50 to +6.00 in 0.50 steps
Material	Methafilcon B 43%, water 57%
Saline Additive	BioMoist
Packaging	30 pack of lenses
Trial lenses available.	



MONTHLY LENSES
3 Pack

Base curve	8.60
Diameter	14.30
Powers available	-0.50 to -6.00 in 0.25 steps -6.00 to -15.00 in 0.50 steps +0.50 to +6.00 in 0.50 steps
Material	Methafilcon A 45%, water 55%
Saline Additive	BioMoist
Packaging	3 pack of lenses
Trial lenses available.	

T 1800 335 559 (Aus only)
E orders@gelflex.com
www.gelflex.com



Gelflex 
Australian contact lenses