

Sofclear ENHANCE

-Fashion Lenses-

— with *BioMoist* —

Model is wearing Aqua #DT105.

The world's only disposable contact lens that truly brightens and enhances the natural colour of your eyes.

Sofclear Enhance uses a carefully formulated combination of special tints to produce the right colour intensity every time. Available in three striking colours, Sofclear Enhance will have everyone talking.

Comfort - it's not just about the look. With all Sofclear contact lenses, comfort is everything. With radial edge technology providing the thinnest of edges, you will go a long way to find a more comfortable contact lens anywhere in the world. Sofclear Enhance lenses are packed in BioMoist.

The BioMoist formula is a hydrating solution containing hydroxypropyl methylcellulose (HPMC) and polyethylene glycol (PEG), which together mimics the natural tear film of the eye. The combination of these superior wetting agents provides a soft, silky feel which helps to optimise comfort, enhancing the overall wearing experience.

Base curve	8.60
Diameter	14.30
Powers available	Plano
	-0.50 to -6.00 in 0.25 steps
	-6.50 to -8.00 in 0.50 steps +0.50 to +6.00 in 0.50 steps
Material	Methafilcon A 45%, water 55%
Saline Additive	BioMoist
Packaging	2 pack of lenses (Monthly)

Trial lenses available in Plano only.



T 1800 335 559 (Aus only)
E orders@gelflex.com
www.gelflex.com



Gelflex 
Australian contact lenses

F1058-04 The content is correct at the time of distribution, September 2021.