

GELFLEX CONTACT LENSES | Our range of products.

Sofclear
COMFORT
with *BioMoist*

Gelflex 
Australian contact lenses



Sofclear

COMFORT with *BioMoist*



Improved initial comfort.



All day moisture retention.



Increased resistance to protein & lipid build up.

EXCEPTIONAL COMFORT ALL DAY LONG

By combining the benefits of Radial Edge technology with the exclusive properties of BioMoist Gelflex has designed a contact lens that provides superior all day comfort.

The BioMoist formula is a hydrating solution containing hydroxypropyl methylcellulose (HPMC) and polyethylene glycol (PEG), which together mimics the natural tear film of the eye. The combination of these superior wetting agents provides a soft, silky feel which helps to optimise comfort upon initial insertion of the lens. The specially formulated BioMoist solution releases moisture when it is needed most, resulting in superior comfort throughout the entire day.

BioMoist also contains unique cleaning agents that help to resist the build-up of protein and lipids on the lens surface, further enhancing the overall wearing experience.

The addition of BioMoist to our comfort range of Sofclear lenses demonstrates our commitment to producing an exceptionally comfortable contact lens.

Sofclear Comfort with BioMoist - it's about *looking good, feeling great* all day long.



Try BioMoist today and
feel the difference!



	MONTHLY	ONE DAY
Base curve	8.60	8.60
Diameter	14.30	14.30
Powers available	-0.50 to -6.00 in 0.25 steps -6.50 to -15.00 in 0.50 steps +0.50 to +6.00 in 0.50 steps	-0.50 to -6.00 in 0.25 steps -6.50 to -15.00 in 0.50 steps +0.50 to +6.00 in 0.50 steps
Material	Methafilcon A 45%, water 55%	Methafilcon B 43%, water 57%
Additive	BioMoist	BioMoist
Packaging	3 or 6 pack of lenses	12 or 30 pack of lenses



ALL DAY COMFORT



MOISTURE RETENTION



RESISTANCE TO BUILD UP

CONTACT LENSES
Australian made and inspired.

Gelflex 
Australian contact lenses

T +61 8 9443 4944 or
1800 998 071 (Australia only)

F +61 8 9443 4147

E info@gelflex.com

52 Mulgul Road, Malaga
Western Australia 6090

www.gelflex.com

CE2460

The content of this booklet is correct at the time of print, July 2020.
F1018-04