

Sofclear COLOURS CLASSIC



Whatever the statement,
Sofclear Colours contact lenses are
the ultimate fashion accessory.

Sofclear Colours are a unique and vibrant range of fashion colour contact lenses.

The colours and designs are inspired by Australia's picturesque landscapes, stunning coastlines and rich sunset hues. Unique patterns are combined to create striking designs that can be worn to enhance any outfit or to reflect a personality.

All lenses are created using state-of-the-art moulding equipment to produce high definition optics.

Featuring our advanced Radial Edge technology, Sofclear Colours contact lenses are designed not only to be seen and admired, but to provide all day comfort to the wearer.



Base curve	8.60
Diameter	14.20
Powers available	Plano -0.50 to -6.00 in 0.25 steps -6.50 to -8.00 in 0.50 steps
Material	Polyhema 57%, water 43%
Packaging	2 pack of lenses (Monthly)

Trial lenses available in Plano only.

T 1800 335 559 (Aus only)
E orders@gelflex.com
www.gelflex.com



Gelflex 
Australian contact lenses